



Forever Medical

## Portable Dental X Ray equipment

### YJ-DXP01



#### Features:

- 1.High Frequency Advances Technology
- 2.Toshiba D-0711 DC Tube.
- 3.Manuual control,Remote control is optional..
- 4.Digitizatin available.
- 5.Computer progran inside
- 6.Power on self test.
- 7.2 year warranty for the whole machine,durable battery,over 200 times for one charging.
- 8.3D lead protection.



Forever Medical

**Parameters:**

Model	YJ-DXP01
Power	400VA
Focal Point	0.4*0.4mm
Voltage	65KV
Current	2mA
High Voltage Supplier	DC 40KHz
Type	I,B
Operation mode	Indirect loading,continuous working
Exposure time	0.1~2.5s

**YJ-DXP02**



**Features:**

- 1.Computor Touch control.
- 2.Toshiba D-045 Tube.
- 3.High Frequency DC.
- 4.Replaceable battery.
- 5.the most smallest dental x ray in China.

**Parameters:**

Model	YJ-DXP02
X Ray Tube	Toshiba D-045
Output Voltage/Current	65kv/2mA
Focus Spot	0.4mm
Anode Angle	16

**TWO models Difference:**

The 01 type is keystroke, original Japanese Toshiba tube, 0.5 lead protection inside the machine, and built-in rechargeable lithium battery.

The 02 type is the LCD touch screen. It has tooth selection, bite chip and digital selection. The original Toshiba ball tube has 0.5 lead protection inside the machine. It is a rechargeable lithium battery, and the 02 type is smaller than the 01 type

## FAQ

Q: What kind of dental product do you trade in ??

A: We are professional supplier of dental x ray unit,dental x ray sensor,dental autoclave,etc.

Q: Your products are original or copy?

A: All products are original made in China, Japan Spare parts.

Q: How can i order?

A: We will make quotation after you send us your purchase plan (including product name, model and quantity). If you agree with the quotation, please send us your company name, address and telephone for goods delivery. We will make proforma invoice and inform you the payment way. Goods delivery details will also be informed accordingly.

Q: What is your payments mode?

A: Mostly payment should be completed by bank transferred T/T 100% in advance.

For every order, we will have to charge 30USD extra to compensate media bank charges. Western Union and Alibaba Trade Assurance are also acceptable. USD, EUR are acceptable.

Q: What is your shipment terms?

A: We often send goods by DHL, FedEx, UPS, TNT, EMS Express; usually 7-15 workdays you can receive goods after send; For the big weight goods, you can choose Air or Sea shipping, usually 30-45 workdays can arrive to your Port.(Exclude unforeseen circumstances, most shipping time is depend on clients country.)

Q: What about the warranty period and after-sale service?

A: Dental equipment, like dental autoclave, x-ray unit, has one year warranty period. If you have problem in using, you can describe the problem details, we will ask technician to give your solution. Free spare parts can be provided if necessary.

Q: Do you supply OEM service?

A: Yes, we can supply OEM service if you have large quantity.